


Trust
Easy Check

Implant Stability Measuring Device

Any Angle

Disposable Vinyl Cap

Similar Osstell Number




Dentium
For Dentists By Dentists

“**You can count on it**”

The value beyond what is seen

CONTACT SYSTEM



Vinyl Cap

Disposable vinyl cap prevents cross infection (optional)



Dentium ISV Scale

	60	65	70	
	Low	Medium	High	
Indication	Implant Risk Monitor ISV	Full Splint IL	Partial	Single
Surgical Protocol		2-Stage	1 or 2-Stage	1-Stage
Restorative Protocol		Traditional	Early	Immediate



Specifications

Product		Easy Check
Main body	Weight(g)	94
	Weight(mm)	26 x 26 x 224
Cradle	Weight(g)	260
	Size(mm)	79 x 111 x 99
Measurement tolerance		Measurement value \pm 3
Measurement angle		0~90°
Measurement time		Less than 3 seconds
Percussive force		Less than 1.0N
Timer		Turns off automatically after one minute
Rated voltage		3.7V / 0.5A
Charger cable		USB C-Type
Battery capacity		420mA
Storage conditions		Temperature 0~40°C / Humidity 30~85%

EYCK- 2021.01[Rev. 0]